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CLSI 評価法：Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard—Second Edition. Clinical and Laboratory Standards Institute, H26-A2 Vol.30 No.14, 2010.

ICSH 評価法：ICSH guidelines for the evaluation of blood cell analysers including those used for differential leucocyte and reticulocyte counting, INTERNATIONAL COUNCIL FOR STANDARDIZATION IN HAEMATOLOGY, WRITING GROUP: C BRIGGS, N CULP, B DAVIS, G.D'ONOFRIO, G ZINI, S.J MACHIN, on behalf of the International Council for Standardization in Hematology (ICSH), International Journal of Laboratory Hematology 36: 613–627, 2014

JSLH 白血球分類参照法：Japanese Society for Laboratory Hematology flow cytometric reference method of determining the differential leukocyte count: external quality assurance using fresh blood sample, Y KAWAI, Y.NAGAI, E.OGAWA, H. KONDO on behalf of the the Standardization Subcommittee for Blood Cell Counting of the Japanese Society for Laboratory Hematology (JSLH), International Journal of Laboratory Hematology 39: 202–222, 2017

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CLSI 白血球分類参照法：Reference Leukocyte (WBC) Differential Count (Proportional) and Evaluation of Instrumental Methods; Approved Standard—Second Edition. Clinical and Laboratory Standards Institute, H20-A2 Vol.27 No.4. 2007

ICSH 白血球分類参照法：Toward a Reference Method for Leukocyte Differential Counts in Blood: Comparison of Three Flow Cytometric Candidate Methods, Mikael Roussel,1* Bruce H. Davis,2 Thierry Fest,1,3,4 Brent L. Wood,5* on behalf of the International Council for Standardization in Hematology (ICSH), Cytometry Part A 81A: 973-982, 2012

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